




AquaVENT® FD140

The AquaVENT FD140 flow driver provides quick and simple transition between high flow oxygen therapy and CPAP

Lifetime paramagnetic O₂ sensor

An interactive LCD touch screen

Integrated drug nebuliser



Intelligent alarm management including high and low pressure, max respiratory rate and apnea alarms

Max Flow rate
CPAP 140L/min,
HFOT 60L/min

60 min internal battery life

Ultra-PEP is used to mobilise secretions, assist with airway clearance and prevent or reverse post-operative atelectasis.

Humidified high flow with **AquaNASE®** maintains mucociliary function and improves oxygenation by delivering flows close to the patients' peak inspiratory requirement.



Ultra-soft nasal prongs

Button clip for extra security

Easy glide neck strap

Armstrong Medical
Wattstown Business Park,
Newbridge Road,
Coleraine, BT52 1BS,
Northern Ireland.

T +44 (0) 28 7035 6029
F +44 (0) 28 7035 6875

info@armstrongmedical.net

armstrongmedical.net



In Stock NOW! Call us for more information

© Armstrong Medical Ltd. 2019